



**Junior Egg Interface**  
Product Code: MN-INT-E

Connects the Junior Egg Logger to a PC or chart reader via serial port



**USB to Serial Port Adapter**  
Product Code: EA-USB-SER

Adapts standard serial port interface to a USB type interface




**Magnetic Key Ring**  
Product Code: EA-KEY

Used when starting the logger as well as time/date stamping to place a bookmark on the reading

## Console Pro Software




- Program and configure logger
- Download & view data analysis
- Export data in CSV & TXT format
- Re-arm function
- Password protected
- MKT, PHI & Growing Degree Day
- CFR 21 part 11



**GRAPH**

Index	Elapsed Time	Date	Time	Internal (°C)
Serial Number				1217-032
H/w version				2.4e
Trip No.				0
Description				
Battery				OK
Start		4/26/2012	4:06:00 PM	
Finish		4/27/2012	8:22:00 AM	
Sampling				1 Minute
Readings				977 (976 Minutes)
Active				-22.2 to 61.3 °C
Specification				-20.0 to 60.0 °C
Average				26.7
Highest				59.8
Lowest				-18.8
PHI Aerobic				5.26
PHI Anaerobic				3.03
MKT				48.60
Growing Degree Day				11.4
Out Of Specification				None
1	00:00:00	4/26/2012	4:06:00 PM	23.2
2	00:01:00	4/26/2012	4:07:00 PM	23.5
3	00:02:00	4/26/2012	4:08:00 PM	23.5
4	00:03:00	4/26/2012	4:09:00 PM	23.7
5	00:04:00	4/26/2012	4:10:00 PM	23.7
6	00:05:00	4/26/2012	4:11:00 PM	23.7
7	00:06:00	4/26/2012	4:12:00 PM	23.7
8	00:07:00	4/26/2012	4:13:00 PM	23.7
9	00:08:00	4/26/2012	4:14:00 PM	23.7

**Readings & Summary**



For more information call us at 732-346-9200 x131, x150 or visit us at [cryopak.com](http://cryopak.com)